

LIFEPAK CR PLUS DEFIBRILLATOR
USER'S CHECKLIST

Instruction	Recommended Corrective Action	Date							
		Initials							
1. Check readiness display for: OK Indicator CHARGE-PAK indicator ATTENTION indicator WRENCH indicator	None. Replace CHARGE-PAK and QUIK-PAK electrode packet. Refer to operating instructions. Contact authorized service personnel.								
2. Check Use By date on QUIK-PAK electrode packet	Replace electrode packet and CHARGE-PAK if date passed.								
3. Check additional supplies.	Replenish as needed.								
4. Check defibrillator for: Damage or cracks Foreign substances	Contact authorized service personnel. Clean the device.								
5. Other:									